VIVADENT®
Vivadent
Technicel

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0132

INSTRUCTIONS FOR USE

MODELLATION
Configure la struttura in composito in forma anatomică ridotta tenendo presente il tipo di rivestimento previsto. La Wandstärke bei Einzelkronen soll mindestens 0,3 mm, aber abutment crowns thickness a minimum of 0.5 mm. Make sure the framework demonstrates adequate stability of shape. Avoid sharp angles. Design the connector areas to be adequate for the position and alloy being used.

SPRUNG
Provide the modeled bridge framework or coping with sprues of a suitable size. Use the direct or in-direct technique being sure that the reservoir is positioned in the heat center. The connection sprue between the reservoir and the coping should be 2.5-3.0 mm in length and width.

INVESTMENT
Weigh the wax pattern including the sprue to determine the quantity of the alloy to be used. (See wax conversion sheet/formula: weight x density = gr. of alloy). Use investment following the manufacturer's instructions.

BURN-OUT
The suggested burnout temperature: 800-850°C/1470-1560°F

MELTING AND CASTING
Use a separate ceramic crucible for each alloy and preheat the crucible in the burnout furnace. Always use fresh alloy. Depending on the type of casting machine, follow the manufacturers instructions for use. If you use Williams Magic Wand, set propane to 0.35 bar/5 psi and oxygen to 0.7 bar/10 psi pressure. Keep the alloy in the reducing atmosphere of the flame between the inner and outer cones. Do not use flux. After casting, bring cool to room temperature. Torch and induction casting: The ingots slumps down, look for oxide break and cast within next 2-4 seconds. Casting Temperature: 1450-1500°C/2642-2732°F

METAL PREPARATION
Carefully clean and clean the object with Al₂O₃. Do not use a hammer for divesting the object to prevent deformation. Finish the framework with carbide burs or with ceramic-bonded grinding instruments. Blast surface with 50-100 micron Al₂O₃ at max. 4.5 bar/65 psi pressure. Subsequently, clean with ultrasonic or steam. Avoid inhalation of dust during grinding!

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SOLDERS AND FLUXES
Design the soldering patry as small as possible and preheat it in the furnace at approximately 600°C/1112°F. The soldering gap should be the same thickness as the soldering strip. Allow the object to cool slowly after soldering.

PRE-Solder: Special High Fusing White Ceramic Solder, High Fusing Bondal Flux NP

Post-Solder: Low Fusing White Gold Solder, High Fusing Bondal Flux NP

Laser weld material: Laser Ceramic White

POLISHING
After glazing remove oxide and flux residue and finish/polish the framework with rubber finishers/polishers.

INDICATIONS
Recommended for crowns, telescope crowns, conus crowns, posts, short- and long span bridges, PMF crowns, implant superstructures, partial dentures.

CONTRAINDICATIONS
For patients with known allergy/sensitivity to any major or minor constituents of this alloy, consultation with a physician is recommended.

SIDE EFFECTS
In individual cases, sensitivity or allergies to components of this alloy may occur.

INTERACTIONS
Galvanic effects may occur between different alloys in the same oral environment.

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SPRUING METHOD

RECOMMENDATIONS
DIRECT: single crowns, inlays and onlays
INDIRECT: multiple units and multiple single crowns

INSTRUCTIONS:
1. Prepare a sprue with a reservoir equal to or larger than the thickest cross-section of the restoration.
2. Maintain the reservoir(s) within the heat center of the investment; keep the restoration(s) approximately 5 mm from the end of the investment and no closer than 5 mm from the sides.
3. Make sure that the sprue is connected to the thickest part of the restoration.
4. The connection between the sprue and the restoration should be flared ("trumpet" configuration); eliminating turbulence of the alloy (causing investment erosion) and facilitating an undisturbed flow of alloy during casting and solidification.
5. Maintain proper alloy weight in order to eliminate the negative effect of a too large button during solidification. The rule of thumb for proper alloy weight determination is:
total wax weight x specific gravity of the alloy = proper alloy weight.

SUGGESTIONS:
1. Thickness and design of wax pattern: follow your manufacturers instructions.
2. Use chill-vents when casting heavy and/or large restorations.

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