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(EN)

GC Initial CAST NP Cobalt-Chrome (CoCr) based casting alloy for metal ceramic veneering (Type 4)

GC Initial CAST NP is a cobalt chrome based casting alloy for fabricating crown & bridges with a suitable metal veneering ceramic or acrylic veneering material. Thanks to its thermal expansion of 14,1 x 10⁻⁶ K⁻¹ it is suitable for the most common metal veneering ceramics on the market. It has outstanding properties like an optimal casting performance, an extra fine molecular structure resulting in smooth, compact surfaces offering an easy and straightforward processing. The alloy offers a reduced oxide formation even after repeated ceramic firings.

GC Initial CAST NP is produced conform the **ISO9693-1:2012** and ISO 22674:2006 standards. It has passed the cytotoxicity test according to the norm ISO 10993-5. The alloy is highly corrosion resistant and it is free of beryllium, indium and gallium.

COMPOSITION

Cobalt	Co	63%
Chromium	Cr	28%
Tungsten	W	3%
Niobium	Nb	4%
Silicon	Si	1%
Others	Fe / Mo	1%

ALLOY CHARACTERISTICS

Solidus-liquidus temperature: 1253 - 1304°C
Casting temperature: 1450 °C
Highest firing temperature: 950°C
Thermal expansion coefficient: (25 - 100°C) 14,1 x 10⁻⁶ K⁻¹ / (25 - 600°C) 14,5 x 10⁻⁶ K⁻¹
Density: 8,4 g/cm³
Vickers Hardness: 324 HV10
Crossbar elongation at fracture: 3,4%
Yield load strength: 475 MPa
Modulus of elasticity: 194 GPa
Colour: White
Recommended ceramics: GC Initial MC

INSTRUCTIONS FOR USE I. WAX-UP

Wax walls should be 0,5 mm thick to ensure that after the manufacturing process the metal wall thickness is at least 0,3 mm. Avoid sharp edges and undercut. Frameworks need to be anatomically reduced ensuring a uniform thickness of the ceramic layer. Connectors should be modelled as strong and as high as possible.

Single crowns: can be directly sprued. Suggested round casting sprues: Ø 3,0 mm to 4,0 mm. Bridge spruing technique is advised. Connections to crossbar Ø 3 to 3,5mm. Distance from modelation to crossbar 2,5mm. Thickness of the crossbars Ø 5mm. Thickness of sprues to crossbar Ø 3 to 4mm.

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II. INVESTING

Use only phosphate bonded investment material for Crown & Bridges. GC Fijvest Premium and GC Fijvest II are optimized for producing high precision fit castings in combination with smooth cast surfaces. Guidelines for optimal use of GC's phosphate bonded investments for Crown & Bridge can be found in the "GC Focus Edition – Phosphate Bonded Investments for Crown & Bridge" and the temperature of 850-950°C is advised. Follow the manufacturer's Instruction for Use related to the pre-heating times before casting.

III. CASTING

We recommend to only use fresh alloys in own crucibles. Sprues and cones can only be re-used once, given that new metal coming from the same lot is added in equal quantity.

Melt GC Initial CAST NP in pre-heated crucibles used for this alloy only.

Amount of casting alloy to use = wax framework weight with sprues / wax density x metal density

High-frequency melting equipment

When last ingots have melted in the pre-heated crucible and the casting immediately start the casting process. Do not use any flux.

Open-flame melting

Place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux.

Overheating the melted material could cause shrinking hollows, micro-porosity and building of a coarse-grained structure and cause bridges to break or ceramic veneering to crack.

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IV. PREPARATION FOR VENEERING

Use tungsten burs to grind and finish the prepared frameworks. Use always the same grinding instrument to avoid contamination. Do not use diamond burs. Grind the surface with slow speed and low pressure. Grind the surface in one direction only to avoid overlaying layers (to avoid air bubbles after the first Opaque firing). Sandblast the prepared framework using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a pressure of 2,5-3 bar., then steam clean. The framework should have a homogeneous grey surface. Proceed with firing process using a compatible veneering ceramic following the manufacturer's instructions.

VI. VENEERING

- Veneering with GC Initial MC: please refer to GC Initial MC Technical Manual
- Veneering with other ceramics: please refer to the manufacturer's instruction for use
- Veneering with acrylics: please refer to the manufacturer's instruction for use.

PACKAGING

GC Initial® CAST NP, 250g
GC Initial® CAST NP, 500g
GC Initial® CAST NP, 1000g

CAUTION

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Before prosthesis application verify if other metal implants are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- Metal dusts and smoke are dangerous for health. Use exhaust fans while casting and sandblasting.
- We recommend to keep batch records and link it to patient's file to allowing its complete traceability.
- This product does not need any special preservation precautions. Nonetheless it has to be employed within the expiry date on the packaging.
- Dispose of all waste according to local regulations.

GUARANTEE

All recommendations are based on GC findings and GC internal tests and are therefore only to be used as guidelines. Any claim for compensation GC receives, in any event, shall not exceed the value of the individual order of GC Initial CAST NP.

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GC Initial CAST NP Cobalt-Chrome (CoCr) based casting alloy for metal ceramic veneering (type 4)

GC Initial CAST NP is a cobalt chrome based casting alloy for fabricating crown & bridges with a suitable metal veneering ceramic or acrylic veneering material. Thanks to its thermal expansion of 14,1 x 10⁻⁶ K⁻¹ it is suitable for the most common metal veneering ceramics on the market. It has outstanding properties like an optimal casting performance, an extra fine microleakage structure resulting in smooth, compact surfaces offering an easy and straightforward processing. The alloy offers a reduced oxide formation even after repeated ceramic firings.

GC Initial CAST NP is produced conform the **ISO 9693-1:2012** and ISO 22674:2006 standards. It has passed the cytotoxicity test according to the norm ISO 10993-5. The alloy is highly corrosion resistant and it is free of beryllium, indium and gallium.

COMPOSITION

Cobalt	Co	63%
Chromium	Cr	28%
Tungsten	W	3%
Niobium	Nb	4%
Silicon	Si	1%
Others	Fe / Mo	1%

ALLOY CHARACTERISTICS
Solidus-liquidus temperature: 1253 - 1304°C
Casting temperature: 1450 °C
Highest firing temperature: 950°C
Modulus of elasticity: 194 GPa
Yield force strength: 475 MPa
Density: 8,4 g/cm³
Vickers hardness: 324 HV10
Elongation at fracture: 3,4%
Tensile force strength: 475 MPa
Modulus of elasticity: 194 GPa
Colour: White
Recommended ceramics: GC Initial MC

INSTRUCTIONS FOR USE I. WAX-UP
Wax walls should be 0,5 mm thick to ensure strength during the manufacturing process the metal wall thickness is at least 0,3 mm. Avoid sharp edges and undercut.
Frameworks need to be anatomically correct with a uniform thickness of the ceramic layer. Connectors should be modelled as strong and as high as possible.

Single crowns:
can be directly sprued. Suggested round casting sprues: Ø 3,0 mm to 4,0 mm.
Bridges:
a crossbar spruing technique is advised. Crossbars to crossbar Ø 3 to 3,5 mm. Distance from modulation to crossbar 2,5mm. Thickness of the crossbars Ø 5mm. Thickness of sprues to crossbar Ø 3 to 4mm.

The higher the volume of the casting the thicker the sprue should be. The casting sprue depends on the object dimensions.

II. INVESTING
Use only phosphate bonded investment material for Crown & Bridges.
GC Initial CAST NP is preadvantitely optimized for producing high precision fit castings in combination with smooth cast surfaces. Guidelines for optimal use of GC's investment material for Crown & Bridge can be found in the "GC Focus Edition – Phosphate Bonded Investments for C&B techniques". A pre-heating temperature of 850-950°C is advised.

Follow the manufacturer's Instruction for Use related to the pre-heating times before casting.
III. CASTING
We recommend to only use fresh alloys in our own crucibles. Sprues and cones can only be re-used once, given that new metal comes from the same lot is added in equal quantity.
Melt GC Initial CAST NP in pre-heated crucibles used for this alloy only.

Amount of casting alloy to use = wax framework weight with sprues / wax density x metal density.
High-frequency melting equipment
When last ingots have melted in the pre-heated ceramic crucible and the metal is fully melted, immediately start the casting process. Do not use any flux.

Open-flame melting
Place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots are melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux.

Overheat the melted material could cause shrinking hollows, micro-porosity and cracks.
Avoid casting in a vacuum, as it may cause bridges to break or ceramic veneering to crack.

Follow the Instructions of the respective manufacturer of the casting devices for parameters and casting procedures.

After casting bench-cool the casting ring. Remove the investment material with care and sandblast using a non-recycling sandblaster with 100 to 125 µm aluminium oxide at a pressure of 3-4 bar.

IV. PREPARATION FOR VENEERING
Use tungsten burs to grind and finish the prepared frameworks. Use always the same cutting and finishing layers to avoid contamination. Do not use diamond burs. Grind the surface with slow speed and low pressure. Grind the surface in one direction only to avoid air bubbles after the first Opaque firing.
Sandblast the prepared framework using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a maximum pressure of 3-4 bar and then steam clean.

V. OXIDE FIRING
Veneer the framework in vacuum at 950-980°C (10°C more than opaque material firing temperature). After firing, carefully sandblast the oxide layer using a non-recycling sandblaster with 110 to 150 µm aluminium oxide at a pressure of 2,5-3 bar, then steam clean. The framework should have a homogeneous grey surface. Proceed with firing process using a compatible veneering ceramic following the manufacturer's instruction.

VI. VENEERING
a) Veneering with GC Initial MC: please refer to GC Initial MC Technical Manual
b) Veneering with other ceramics: please refer to the manufacturer's instruction for use
c) Veneering with acrylics: please refer to the manufacturer's instruction for use.

PACKAGING
GC Initial[®] CAST NP, 250g
GC Initial[®] CAST NP, 500g
GC Initial[®] CAST NP, 1000g

CAUTION
1) Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects.
2) Before prostheses application verify if other metal implants are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
3) Metal dust and smoke are dangerous for health. Use exhaust fans while casting and sandblasting.
4) We recommend to keep batch records for each patient's file to allowing its complete traceability.
5) This product does not need any special preservation precautions. Nonetheless it is recommended to observe the expiry date on the packaging.
6) Dispose of all waste according to local regulations.

GUARANTEE
All recommendations are based on GC findings and GC internal tests and are subject to any change in guidelines. Any claim for compensation GC receives, in any event, shall not exceed the value of the individual order of GC Initial CAST NP.

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